

For Research Use Only

# CCDC91 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84010-3-PBS

## Basic Information

Catalog Number: 84010-3-PBS	GenBank Accession Number: BC028682	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55297	CloneNo.: 241159E2
Source: Rabbit	UNIPROT ID: Q7Z6B0	
Isotype: IgG	Full Name: coiled-coil domain containing 91	
Immunogen Catalog Number: AG34863		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00946-2, CCDC91 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84010-2-PBS. Detection antibody: 84010-3-PBS. Standard: Ag34863. Range: 0.313-40 ng/mL.