

For Research Use Only

ABR Recombinant antibody, PBS Only (Detector)



Catalog Number: 84011-2-PBS

Basic Information

Catalog Number: 84011-2-PBS	GenBank Accession Number: NM_021962	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 29	CloneNo.: 241121F10
Source: Rabbit	UNIPROT ID: Q12979	
Isotype: IgG	Full Name: active BCR-related gene	
Immunogen Catalog Number: AG29522	Calculated MW: 98 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00961-3, ABR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84011-1-PBS. Detection antibody: 84011-2-PBS. Standard: Ag29522. Range: 0.313-40 ng/mL