

For Research Use Only

SRP54 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84014-4-PBS



Basic Information

Catalog Number: 84014-4-PBS	GenBank Accession Number: BC003389	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6729	CloneNo.: 241125H7
Source: Rabbit	UNIPROT ID: P61011	
Isotype: IgG	Full Name: signal recognition particle 54kDa	
Immunogen Catalog Number: AG2327	Calculated MW: 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00955-1, SRP54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84014-4-PBS. Detection antibody: 84014-2-PBS. Standard: Ag2327. Range: 0.78-100 ng/mL