

For Research Use Only

IGF1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 84015-2-PBS

Basic Information

Catalog Number: 84015-2-PBS	GenBank Accession Number: NM_001111285.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3479	CloneNo.: 241094G7
Source: Rabbit	UNIPROT ID: P05019-1	
Isotype: IgG	Full Name: insulin-like growth factor 1 (somatomedin C)	
Immunogen Catalog Number: AG29197	Calculated MW: 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00956-1, IGF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84015-2-PBS. Detection antibody: 84015-3-PBS. Standard: Ag29197. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00956-2, IGF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84015-1-PBS. Detection antibody: 84015-2-PBS. Standard: Ag29197. Range: 0.781-100 ng/mL.