

For Research Use Only

SUMO3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84017-2-PBS



Basic Information

Catalog Number: 84017-2-PBS	GenBank Accession Number: BC008420	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6612	CloneNo.: 241156C3
Source: Rabbit	UNIPROT ID: P55854	
Isotype: IgG	Full Name: SMT3 suppressor of mif two 3 homolog 3 (S. cerevisiae)	
Immunogen Catalog Number: AG1135	Calculated MW: 12 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00968-1, SUMO2/3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84017-2-PBS. Detection antibody: 84017-1-PBS. Standard: Ag1135. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84017-2-PBS against Human SUMO3 was performed. The affinity constant is 2.04 nM.