

For Research Use Only

SUMO3 Recombinant antibody, PBS Only (Capture)

Catalog Number:84017-6-PBS



Basic Information

Catalog Number: 84017-6-PBS	GenBank Accession Number: BC008420	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6612	CloneNo.: 241156E12
Source: Rabbit	UNIPROT ID: P55854	
Isotype: IgG	Full Name: SMT3 suppressor of mif two 3 homolog 3 (S. cerevisiae)	
Immunogen Catalog Number: AG1135	Calculated MW: 12 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00968-3, Human SUMO3 Recombinant Matched Antibody Pair - PBS only. 84017-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1135. 84017-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL

The mean SUMO3 concentration was determined to be 47.49 ng/mL in K-562 cell extract based on a 2.4 mg/mL extract load.