

For Research Use Only

FAP Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 84018-3-PBS

Basic Information

Catalog Number: 84018-3-PBS	GenBank Accession Number: NM_004460.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2191	CloneNo.: 241129F6
Source: Rabbit	UNIPROT ID: Q12884-1	
Isotype: IgG	Full Name: fibroblast activation protein, alpha	
	Calculated MW: 88 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00971-1, FAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84018-2-PBS. Detection antibody: 84018-3-PBS. Standard: Eg0821. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00971-3, FAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84018-3-PBS. Detection antibody: 84018-1-PBS. Standard: Eg0821. Range: 0.78-50 ng/mL