

For Research Use Only

FOXF1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84020-4-PBS



Basic Information

Catalog Number: 84020-4-PBS	GenBank Accession Number: NP_001442	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 2294	CloneNo.: 241076D1
Source: Rabbit	UNIPROT ID: Q12946	
Isotype: IgG	Full Name: forkhead box F1	
Immunogen Catalog Number: AG34266	Calculated MW: 40 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00943-2, Human FOXF1 Recombinant Matched Antibody Pair - PBS only. 84020-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34266. 84020-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.



The mean FOXF1 concentration was determined to be 9.82 ng/mL in HCT 116 cell extract based on a 4.80 mg/mL extract load.