

For Research Use Only

# NFKB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84022-3-PBS

## Basic Information

<b>Catalog Number:</b> 84022-3-PBS	<b>GenBank Accession Number:</b> BC002844	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4791	<b>CloneNo.:</b> 241126F6
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000077150	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q00653	
<b>Immunogen Catalog Number:</b> AG0662	<b>Full Name:</b> nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)	
	<b>Calculated MW:</b> 97 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00950-2, NFkB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84022-3-PBS. Detection antibody: 84022-4-PBS. Standard: Ag0662. Range: 0.781-100 ng/mL.