

For Research Use Only

CSF2RA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84023-3-PBS



Basic Information

Catalog Number:

84023-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC002635

GeneID (NCBI):

1438

UNIPROT ID:

P15509

Full Name:

colony stimulating factor 2 receptor, alpha, low-affinity (granulocyte-macrophage)

Calculated MW:

46 kDa

Purification Method:

Protein A purification

CloneNo.:

241266C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00953-2, CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84023-4-PBS. Detection antibody: 84023-3-PBS. Standard: RP02335. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00953-3, CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84023-2-PBS. Detection antibody: 84023-3-PBS. Standard: RP02335. Range: 0.625-80 ng/mL.