

For Research Use Only

TLR4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84025-1-PBS

Basic Information

Catalog Number: 84025-1-PBS	GenBank Accession Number: NM_138554	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7099	CloneNo.: 241272A11
Source: Rabbit	UNIPROT ID: O00206	
Isotype: IgG	Full Name: toll-like receptor 4	
	Calculated MW: 96 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00959-1, TLR4/CD284 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84025-2-PBS. Detection antibody: 84025-1-PBS. Standard: RP02344. Range: 6.25-200 ng/mL.