

For Research Use Only

TLR4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84025-2-PBS

Basic Information

Catalog Number:

84025-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_138554

GeneID (NCBI):

7099

UNIPROT ID:

O00206

Full Name:

toll-like receptor 4

Calculated MW:

96 kDa

Purification Method:

Protein A purification

CloneNo.:

241272E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00959-1, TLR4/CD284 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84025-2-PBS. Detection antibody: 84025-1-PBS. Standard: RP02344. Range: 6.25-200 ng/mL.