

For Research Use Only

ACOX3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84027-1-PBS

Basic Information

Catalog Number: 84027-1-PBS	GenBank Accession Number: BC017053	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8310	CloneNo.: 241172A11
Source: Rabbit	UNIPROT ID: O15254	
Isotype: IgG	Full Name: acyl-Coenzyme A oxidase 3, pristanoyl	
Immunogen Catalog Number: AG10553	Calculated MW: 624 aa, 70 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00962-1, ACOX3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84027-3-PBS. Detection antibody: 84027-1-PBS. Standard: Ag10553. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00962-2, ACOX3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84027-2-PBS. Detection antibody: 84027-1-PBS. Standard: Ag10553. Range: 1.563-100 ng/mL.