

For Research Use Only

TMPPE Recombinant antibody, PBS Only (Detector)



Catalog Number: 84028-1-PBS

Basic Information

Catalog Number: 84028-1-PBS	GenBank Accession Number: BC136873	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 643853	CloneNo.: 241180B10
Source: Rabbit	UNIPROT ID: Q6ZT21	
Isotype: IgG	Full Name: transmembrane protein with metallophosphoesterase domain	
Immunogen Catalog Number: AG19556	Calculated MW: 453 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00963-1, TMPPE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84028-2-PBS. Detection antibody: 84028-1-PBS. Standard: Ag19556. Range: 0.313-40 ng/mL