

For Research Use Only

# PEG3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84030-1-PBS

## Basic Information

Catalog Number: 84030-1-PBS	GenBank Accession Number: BC150272	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5178	CloneNo.: 24123485
Source: Rabbit	UNIPROT ID: Q9GZU2	
Isotype: IgG	Full Name: paternally expressed 3	
Immunogen Catalog Number: AG20132	Calculated MW: 1588 aa, 181 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00940-3, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-1-PBS. Detection antibody: 84030-4-PBS. Standard: Ag20132. Range: 0.156-20 ng/mL