

For Research Use Only

PEG3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84030-4-PBS

Basic Information

Catalog Number: 84030-4-PBS	GenBank Accession Number: BC150272	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5178	CloneNo.: 241234G1
Source: Rabbit	UNIPROT ID: Q9GZU2	
Isotype: IgG	Full Name: paternally expressed 3	
Immunogen Catalog Number: AG20132	Calculated MW: 1588 aa, 181 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.com

W: ptgcn.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP00940-2, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-2-PBS. Detection antibody: 84030-4-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP00940-2, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-2-PBS. Detection antibody: 84030-4-PBS. Standard: Ag20132. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP00940-3, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-1-PBS. Detection antibody: 84030-4-PBS. Standard: Ag20132. Range: 0.156-20 ng/mL.