

For Research Use Only

ACSF2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84033-1-PBS

Basic Information

Catalog Number: 84033-1-PBS	GenBank Accession Number: BC012053	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 80221	CloneNo.: 241208C1
Source: Rabbit	UNIPROT ID: Q96CM8	
Isotype: IgG	Full Name: acyl-CoA synthetase family member	
Immunogen Catalog Number: AG9199	Calculated MW: 615 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00952-2, ACSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84033-1-PBS. Detection antibody: 84033-2-PBS. Standard: Ag9199. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00952-3, ACSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84033-1-PBS. Detection antibody: 84033-3-PBS. Standard: Ag9199. Range: 2.5-80 ng/mL