

For Research Use Only

ACSF2 Recombinant antibody

Catalog Number: 84033-3-RR



Basic Information

Catalog Number:

84033-3-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9199

GenBank Accession Number:

BC012053

GeneID (NCBI):

80221

UNIPROT ID:

Q96CM8

Full Name:

acyl-CoA synthetase family member

2

Calculated MW:

615 aa, 68 kDa

Observed MW:

51-71 kDa

Purification Method:

Protein A purification

CloneNo.:

241208D9

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, SGC-7901 cells

Background Information

ACSF2 (encoded by acyl-CoA synthetase family member 2) belongs to the acyl-CoA synthetase (ACS) family, activating fatty acids by forming a thioester bond with CoA.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84033-3-RR (ACSF2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.