

For Research Use Only

RIPK1-Specific Recombinant antibody, PBS Only (Detector)



Catalog Number: 84036-1-PBS

Basic Information

Catalog Number: 84036-1-PBS	GenBank Accession Number: NM_003804	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8737	CloneNo.: 241099C7
Source: Rabbit	UNIPROT ID: Q13546	
Isotype: IgG	Full Name: receptor (TNFRSF)-interacting serine-threonine kinase 1	
	Calculated MW: 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00944-1, RIPK1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84036-2-PBS. Detection antibody: 84036-1-PBS. Standard: Ag31943. Range: 1.563-100 ng/mL