

For Research Use Only

Mouse CD98/SLC3A2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84040-1-PBS



Basic Information

Catalog Number:

84040-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008577.4

GeneID (NCBI):

17254

UNIPROT ID:

P10852-1

Full Name:

solute carrier family 3 (activators of dibasic and neutral amino acid transport), member 2

Calculated MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

241247A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00969-2, MOUSE CD98/SLC3A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84040-1-PBS. Detection antibody: 84040-3-PBS. Standard: Eg1374. Range: 1.25-40 ng/mL



Cytometric bead array standard curve of MP00969-3, MOUSE CD98/SLC3A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84040-1-PBS. Detection antibody: 84040-5-PBS. Standard: Eg1374. Range: 1.25-40 ng/mL