

For Research Use Only

PRC1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84041-2-PBS

Basic Information

Catalog Number: 84041-2-PBS	GenBank Accession Number: BC003138	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9055	CloneNo.: 241216H2
Source: Rabbit	UNIPROT ID: O43663	
Isotype: IgG	Full Name: protein regulator of cytokinesis 1	
Immunogen Catalog Number: AG8028	Calculated MW: 72 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00970-1, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-1-PBS. Standard: Ag8028. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard: null. Range: null.



Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard: Ag8028. Range: 0.625-80 ng/mL