

For Research Use Only

GPAM Recombinant antibody, PBS Only (Detector)



Catalog Number: 84054-1-PBS

Basic Information

Catalog Number: 84054-1-PBS	GenBank Accession Number: BC030783	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57678	CloneNo.: 241144B3
Source: Rabbit	UNIPROT ID: Q9HCL2	
Isotype: IgG	Full Name: glycerol-3-phosphate acyltransferase, mitochondrial	
Immunogen Catalog Number: AG3120	Calculated MW: 828 aa, 94 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00994-2, GPAM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84054-5-PBS. Detection antibody: 84054-1-PBS. Standard: Ag3120. Range: 0.625-80 ng/mL