

For Research Use Only

MIRO1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84055-4-PBS



Basic Information

Catalog Number: 84055-4-PBS	GenBank Accession Number: BC060781	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55288	CloneNo.: 241157H2
Source: Rabbit	UNIPROT ID: Q8IXI2	
Isotype: IgG	Full Name: ras homolog gene family, member T1	
Immunogen Catalog Number: AG15605	Calculated MW: 618 aa, 71 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01000-3, MIRO1/RHOT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84055-2-PBS. Detection antibody: 84055-4-PBS. Standard: Ag15605. Range: 1.25-40 ng/mL.