

For Research Use Only

RHOT1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84055-7-PBS



Basic Information

Catalog Number: 84055-7-PBS	GenBank Accession Number: BC060781	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55288	CloneNo.: 241157A9
Source: Rabbit	UNIPROT ID: Q8IXI2	
Isotype: IgG	Full Name: ras homolog gene family, member T1	
Immunogen Catalog Number: AG15605	Calculated MW: 618 aa, 71 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01000-4, Human MIRO1 Recombinant Matched Antibody Pair - PBS only. 84055-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag15605. 84055-6-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL.



The mean MIRO1 concentration was determined to be 17.32 ng/mL in mouse brain cell extract based on a 4.60 mg/mL extract load and 25.35 ng/mL in HeLa cell extract based on a 4.60 mg/mL extract load.