

For Research Use Only

CASQ2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84056-1-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 84056-1-PBS | GenBank Accession Number: BC022288 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 845 | CloneNo.: 241211B2 |
| Source: Rabbit | UNIPROT ID: O14958 | |
| Isotype: IgG | Full Name: calsequestrin 2 (cardiac muscle) | |
| Immunogen Catalog Number: AG33584 | Calculated MW: 46 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01004-2, CASQ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84056-5-PBS. Detection antibody: 84056-1-PBS. Standard: Ag33584. Range: 0.313-40 ng/mL