

For Research Use Only

RNF40 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84058-2-PBS

Basic Information

Catalog Number: 84058-2-PBS	GenBank Accession Number: BC006133	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9810	CloneNo.: 241214C8
Source: Rabbit	UNIPROT ID: O75150	
Isotype: IgG	Full Name: ring finger protein 40	
Immunogen Catalog Number: AG8047	Calculated MW: 114 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00974-2, RNF40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84058-2-PBS. Detection antibody: 84058-3-PBS. Standard: Ag8047. Range: 0.156-20 ng/mL