

For Research Use Only

KLK3/PSA Recombinant antibody, PBS Only (Detector)



Catalog Number: 84059-3-PBS

Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 84059-3-PBS | GenBank Accession Number: NM_001648.2 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 354 | CloneNo.: 241223G9 |
| Source: Rabbit | UNIPROT ID: P07288-1 | |
| Isotype: IgG | Full Name: kallikrein-related peptidase 3 | |
| | Calculated MW: 29 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00975-1, KLK3/PSA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84059-2-PBS. Detection antibody: 84059-3-PBS. Standard: Eg0869. Range: 0.156-20 ng/mL.