

For Research Use Only

SLC9A3R2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84060-1-PBS

Basic Information

Catalog Number: 84060-1-PBS	GenBank Accession Number: BC106001	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9351	CloneNo.: 241215B1
Source: Rabbit	UNIPROT ID: Q15599	
Isotype: IgG	Full Name: solute carrier family 9 (sodium/hydrogen exchanger), member 3 regulator 2	
Immunogen Catalog Number: AG13476	Calculated MW: 337 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00976-1, SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84060-4-PBS. Detection antibody: 84060-1-PBS. Standard: Ag13476. Range: 0.313-40 ng/mL.