

For Research Use Only

SRP19 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84061-1-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 84061-1-PBS | GenBank Accession Number: BC010947 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 6728 | CloneNo.: 241219C12 |
| Source: Rabbit | UNIPROT ID: P09132 | |
| Isotype: IgG | Full Name: signal recognition particle 19kDa | |
| Immunogen Catalog Number: AG8903 | Calculated MW: 144 aa, 16 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00977-1, SRP19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84061-1-PBS. Detection antibody: 84061-2-PBS. Standard: Ag8903. Range: 1.25-40 ng/mL