

For Research Use Only

LPAR5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84065-2-PBS

Basic Information

Catalog Number: 84065-2-PBS	GenBank Accession Number: BC033571	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57121	CloneNo.: 241221F11
Source: Rabbit	UNIPROT ID: Q9H1C0	
Isotype: IgG	Full Name: lysophosphatidic acid receptor 5	
Immunogen Catalog Number: AG24253	Calculated MW: 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00981-1, LPAR5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84065-2-PBS. Detection antibody: 84065-1-PBS. Standard: Ag24253. Range: 0.313-40 ng/mL