

For Research Use Only

Mouse CD138/Syndecan-1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84069-4-PBS



Basic Information

Catalog Number: 84069-4-PBS	GenBank Accession Number: NM_011519.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 20969	CloneNo.: 241278G10
Source: Rabbit	UNIPROT ID: P18828-1	
Isotype: IgG	Full Name: syndecan 1	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00990-3, MOUSE CD138/Syndecan-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84069-2-PBS. Detection antibody: 84069-4-PBS. Standard: Eg1221. Range: 0.625-80 ng/mL