

For Research Use Only

CD38 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84070-4-PBS

Basic Information

Catalog Number: 84070-4-PBS	GenBank Accession Number: BC007964	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 952	CloneNo.: 241246E5
Source: Rabbit	ENSEMBL Gene ID: ENSG00000004468	
Isotype: IgG	UNIPROT ID: P28907	
	Full Name: CD38 molecule	
	Calculated MW: 300 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00991-1, CD38 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84070-4-PBS. Detection antibody: 84070-2-PBS. Standard: Eg0110. Range: 0.313-40 ng/mL