

For Research Use Only

PRG3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84075-1-PBS

Basic Information

Catalog Number: 84075-1-PBS	GenBank Accession Number: BC101613	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10394	CloneNo.: 241176A11
Source: Rabbit	UNIPROT ID: Q9Y2Y8	
Isotype: IgG	Full Name: proteoglycan 3	
Immunogen Catalog Number: AG12706	Calculated MW: 225 aa, 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01017-1, PRG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84075-1-PBS. Detection antibody: 84075-4-PBS. Standard: Ag12706. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01017-2, PRG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84075-1-PBS. Detection antibody: 84075-3-PBS. Standard: Ag12706. Range: 0.313-40 ng/mL