

For Research Use Only

FKBP2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84077-3-PBS

Basic Information

Catalog Number: 84077-3-PBS	GenBank Accession Number: BC003384	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2286	CloneNo.: 241166E4
Source: Rabbit	UNIPROT ID: P26885	
Isotype: IgG	Full Name: FK506 binding protein 2, 13kDa	
Immunogen Catalog Number: AG2255	Calculated MW: 142 aa, 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00984-3, FKBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84077-3-PBS. Detection antibody: 84077-4-PBS. Standard: Ag2255. Range: 0.781-100 ng/mL.