

For Research Use Only

CD36/SR-B3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84078-2-PBS

Basic Information

Catalog Number: 84078-2-PBS	GenBank Accession Number: NM_001001548.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 948	CloneNo.: 241262C7
Source: Rabbit	UNIPROT ID: P16671-1	
Isotype: IgG	Full Name: CD36 molecule (thrombospondin receptor)	
	Calculated MW: 53kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00992-2, CD36/SR-B3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84078-2-PBS. Detection antibody: 84078-3-PBS. Standard: Eg1070. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP00992-3, CD36/SR-B3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84078-2-PBS. Detection antibody: 84078-1-PBS. Standard: Eg1070. Range: 6.25-100 ng/mL.