

For Research Use Only

ROR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84079-2-PBS

Basic Information

Catalog Number: 84079-2-PBS	GenBank Accession Number: NM_005012	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4919	CloneNo.: 241264B7
Source: Rabbit	UNIPROT ID: Q01973	
Isotype: IgG	Full Name: receptor tyrosine kinase-like orphan receptor 1	
	Calculated MW: 104kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00987-1, ROR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84079-2-PBS. Detection antibody: 84079-1-PBS. Standard: Eg0203. Range: 0.625-40 ng/mL