

For Research Use Only

# ROR1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84079-4-PBS



## Basic Information

Catalog Number:

84079-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005012

GeneID (NCBI):

4919

UNIPROT ID:

Q01973

Full Name:

receptor tyrosine kinase-like orphan  
receptor 1

Calculated MW:

104kd

Purification Method:

Protein A purification

CloneNo.:

241264D7

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP00987-2, Human ROR1 Recombinant Matched Antibody Pair - PBS only. 84079-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0203. 84079-4-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL

The mean ROR1 concentration was determined to be 16.38 ng/mL in Caco-2 cell extract based on a 1.5 mg/mL extract load.