

For Research Use Only

SEMA3B Recombinant antibody, PBS Only (Detector)



Catalog Number: 84080-2-PBS

Basic Information

Catalog Number: 84080-2-PBS	GenBank Accession Number: BC009113	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7869	CloneNo.: 241171D2
Source: Rabbit	UNIPROT ID: Q13214	
Isotype: IgG	Full Name: sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3B	
Immunogen Catalog Number: AG25163		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00995-1, SEMA3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84080-1-PBS. Detection antibody: 84080-2-PBS. Standard: Ag25163. Range: 0.625-80 ng/mL.