

For Research Use Only

PCF11 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84081-3-PBS

Basic Information

Catalog Number: 84081-3-PBS	GenBank Accession Number: BC146778	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51585	CloneNo.: 241213E8
Source: Rabbit	UNIPROT ID: O94913	
Isotype: IgG	Full Name: PCF11, cleavage and polyadenylation factor subunit, homolog (S. cerevisiae)	
Immunogen Catalog Number: AG20263	Calculated MW: 1555 aa, 173 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00996-1, PCF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84081-2-PBS. Detection antibody: 84081-3-PBS. Standard: Ag20263. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00996-2, PCF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84081-1-PBS. Detection antibody: 84081-3-PBS. Standard: Ag20263. Range: 0.625-80 ng/mL