

For Research Use Only

TRP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84082-4-PBS



Basic Information

Catalog Number: 84082-4-PBS	GenBank Accession Number: BC052608	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7306	CloneNo.: 241242C2
Source: Rabbit	UNIPROT ID: P17643	
Isotype: IgG	Full Name: tyrosinase-related protein 1	
Immunogen Catalog Number: AG10166	Calculated MW: 537 aa, 61 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00997-3, Human TRP1 Recombinant Matched Antibody Pair - PBS only. 84082-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10166. 84082-2-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL.



The mean TRP1 concentration was determined to be 16.83 ng/mL in A431 cell extract based on a 4.00 mg/mL extract load and 27.93 ng/mL in A375 cell extract based on a 4.60 mg/mL extract load.