

For Research Use Only

MYO6 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 84085-3-PBS

Basic Information

Catalog Number: 84085-3-PBS	GenBank Accession Number: BC146764	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4646	CloneNo.: 241251G11
Source: Rabbit	UNIPROT ID: Q9UM54	
Isotype: IgG	Full Name: myosin VI	
Immunogen Catalog Number: AG24906	Calculated MW: 1285 aa, 149 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01001-1, MYO6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84085-3-PBS. Detection antibody: 84085-1-PBS. Standard: Ag24906. Range: 0.781-100 ng/ml.



Cytometric bead array standard curve of MP01001-3, MYO6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84085-2-PBS. Detection antibody: 84085-3-PBS. Standard: Ag24906. Range: 0.781-100 ng/ml.