

For Research Use Only

TCF1/TCF7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84089-1-PBS



Basic Information

Catalog Number:

84089-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5792

GenBank Accession Number:

BC048769

GeneID (NCBI):

6932

ENSEMBL Gene ID:

ENSG00000081059

UNIPROT ID:

P36402

Full Name:

transcription factor 7 (T-cell specific, HMG-box)

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

24115485

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01006-1, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-1-PBS. Detection antibody: 84089-3-PBS. Standard: Ag5792. Range: 1.56-20 ng/mL.



Cytometric bead array standard curve of MP01006-2, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-1-PBS. Detection antibody: 84089-2-PBS. Standard: Ag5792. Range: 1.56-10 ng/mL.