

For Research Use Only

TCF1/TCF7 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84089-2-PBS



Basic Information

Catalog Number: 84089-2-PBS	GenBank Accession Number: BC048769	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6932	CloneNo.: 241154B11
Source: Rabbit	ENSEMBL Gene ID: ENSG00000081059	
Isotype: IgG	UNIPROT ID: P36402	
Immunogen Catalog Number: AG5792	Full Name: transcription factor 7 (T-cell specific, HMG-box)	
	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01006-2, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-1-PBS. Detection antibody: 84089-2-PBS. Standard: Ag5792. Range: 1.56-10 ng/mL.



Cytometric bead array standard curve of MP01006-3, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-2-PBS. Detection antibody: 84089-3-PBS. Standard: Ag5792. Range: 1.56-20 ng/mL.