

For Research Use Only

MYO9B Recombinant antibody, PBS Only (Capture)

Catalog Number:84090-1-PBS



Basic Information

Catalog Number: 84090-1-PBS	GenBank Accession Number: BC018108	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4650	CloneNo.: 241163A5
Source: Rabbit	UNIPROT ID: Q13459	
Isotype: IgG	Full Name: myosin IXB	
Immunogen Catalog Number: AG33349	Calculated MW: 2023 aa, 223 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01007-1, MYO9B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84090-1-PBS. Detection antibody: 84090-3-PBS. Standard: Ag33349. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01007-2, MYO9B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84090-1-PBS. Detection antibody: 84090-2-PBS. Standard: Ag33349. Range: 3.125-200 ng/mL.