

For Research Use Only

MYO9B Recombinant antibody, PBS Only (Capture)

Catalog Number:84090-5-PBS



Basic Information

Catalog Number: 84090-5-PBS	GenBank Accession Number: BC018108	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4650	CloneNo.: 241163B6
Source: Rabbit	UNIPROT ID: Q13459	
Isotype: IgG	Full Name: myosin IXB	
Immunogen Catalog Number: AG33349	Calculated MW: 2023 aa, 223 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01007-3, Human MYO9B Recombinant Matched Antibody Pair - PBS only. 84090-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag33349. 84090-3-PBS was HRP conjugated as the detection antibody. Range: 1.56-50 ng/mL



The mean MYO9B concentration was determined to be 38.37 ng/mL in K-562 cell extract based on a 2.10 mg/mL extract load.