

For Research Use Only

CD209 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84095-4-PBS

Basic Information

Catalog Number: 84095-4-PBS	GenBank Accession Number: BC110615	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 30835	CloneNo.: 240907D7
Source: Rabbit	UNIPROT ID: Q9NNX6	
Isotype: IgG	Full Name: CD209 molecule	
	Calculated MW: 404 aa, 45.7 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01013-3, CD209/DC-SIGN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84095-5-PBS. Detection antibody: 84095-4-PBS. Standard: Eg0062. Range: 0.313-20 ng/mL.