

For Research Use Only

CD209 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84095-5-PBS

Basic Information

Catalog Number:

84095-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC110615

GeneID (NCBI):

30835

UNIPROT ID:

Q9NNX6

Full Name:

CD209 molecule

Calculated MW:

404 aa, 45.7 kDa

Purification Method:

Protein A purification

CloneNo.:

240907E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01013-3, CD209/DC-SIGN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84095-5-PBS. Detection antibody: 84095-4-PBS. Standard: Eg0062. Range: 0.313-20 ng/mL.