

For Research Use Only

CFB Recombinant antibody, PBS Only (Detector)



Catalog Number: 84101-2-PBS

Basic Information

Catalog Number: 84101-2-PBS	GenBank Accession Number: BC007990	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 629	CloneNo.: 241141G4
Source: Rabbit	UNIPROT ID: P00751	
Isotype: IgG	Full Name: complement factor B	
Immunogen Catalog Number: AG0226	Calculated MW: 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01031-2, Complement factor B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84101-1-PBS. Detection antibody: 84101-2-PBS. Standard: Ag0226. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01031-3, Complement factor B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84101-3-PBS. Detection antibody: 84101-2-PBS. Standard: Ag0226. Range: 0.781-100 ng/mL.