

For Research Use Only

DDX24 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84103-4-PBS

Basic Information

Catalog Number: 84103-4-PBS	GenBank Accession Number: BC008847	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57062	CloneNo.: 241145H7
Source: Rabbit	UNIPROT ID: Q9GZR7	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 24	
Immunogen Catalog Number: AG8407	Calculated MW: 96 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01040-1, DDX24 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84103-1-PBS. Detection antibody: 84103-4-PBS. Standard: Ag8407. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01040-3, DDX24 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84103-2-PBS. Detection antibody: 84103-4-PBS. Standard: Ag8407. Range: 0.625-40 ng/mL