

For Research Use Only

CCDC134 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84106-4-PBS

Basic Information

Catalog Number: 84106-4-PBS	GenBank Accession Number: BC017693	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79879	CloneNo.: 241212H5
Source: Rabbit	UNIPROT ID: Q9H6E4	
Isotype: IgG	Full Name: coiled-coil domain containing 134	
Immunogen Catalog Number: AG21990	Calculated MW: 229 aa, 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01057-2, CCDC134 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84106-4-PBS. Detection antibody: 84106-1-PBS. Standard: Ag21990. Range: 0.625-80 ng/mL.